




Die and Model 2 Gray

Rapidly 3D Print Models for Aligner Production



SRI-0202011

SKU Number

SprintRay Die & Model 2 represents a big step forward in dental modelling materials. Our new gray color offers extremely high color contrast to help illustrate details. This new resin features almost no pigment settling, which means you can spend more time printing and less time mixing. Improved adhesion to the build plate is matched by an increase in flexural modulus and strength, providing higher success rates and a reduced chance of breakage.

A reduced odor improves the quality of life for 3D printing models for an overall experience that your practice staff will love. Die & Model 2 Gray features an improved color tone, making it easier than ever to use digital technology to perform advanced restorative dentistry. With rapid, repeatable results that are easier than ever to use, Die & Model 2 Gray is helping to bring dental 3D printing into dental practices everywhere.

KEY BENEFITS

- ✓ Improved flexural modulus and strength
- ✓ Up to 10% faster print speeds
- ✓ Print up to eight full arch models in just 15 minutes*
- ✓ High-contrast color for maximum visibility of small details
- ✓ Easy to separate from thermoforming materials

FLEXURAL STRENGTH

Higher flexural strength ensures that the model doesn't break during normal use.

Die and Model 2 Gray 91.5 MPa

Die and Model Gray 66.7 MPa

FLEXURAL MODULUS

Higher flexural modulus prevents the model from distorting, which can impact accuracy.

Die and Model 2 Gray 2650 MPa

Die and Model Gray 1700 MPa

* When printed on SprintRay Pro with Ludicrous Speed.